




Portable X-ray machine



Used in hospitals | medical clinics | nursing homes

Introduction

Compact and convenient X-ray

YSX050-G is a multifunctional portable X-ray machine, including a 5kw X-ray power supply device. Small size, light weight and high power density, It can fulfill various portable X-ray application scenarios.




Product features



① Compact and portable, Easy to use



② Wireless controlled exposure



③ Wired manual exposure




④ Flexible operation, fast imaging



⑤ Bluetooth Support
Easy connection to computer



Application



Clinical Images



Technical Parameters

Model	YJ-XP50
System Requirements	
Electric source (V)	AC 220 - 240 / AC 100-130 50Hz/60Hz
Power (kW)	5
Generator	
kV range (kVp)	40 -125
mAs range (mAs)	1-100
mA range (mA)	10 -100
Ranging mode	Tapeline
Communication mode	RS232/Bluetooth
Exposure mode	Wired/Wireless
X-ray Tube	
Focal spots (mm)	0.6/1.8
Anode Angle (°)	15
Frame	
Dimension (mm³)	349 x 208 x 406
Weight (kg)	17.7

Dimension

(Unit: mm)

